

Annex B CSOP Pelvic Trauma - T-POD

T-POD Application Instructions

Important Notes

1. The T-POD Must be applied directly to the skin, and not over clothing.
2. Be careful not to over-tension the T-POD

INSTRUCTIONS AND USE


T-POD[®]
Pelvic Stabilizer




1 Slide Belt under patient and into position under the pelvis.




2 Trim the Belt, leaving a 6-8" gap over the centre of the abdomen.



3 Apply Velcro-backed Pulley System on each side of the gap.



4 Draw the Pull Tab, creating simultaneous circumferential compression.



5 Secure the Velcro-backed Pull Tab to the Belt.



6 Record the date and time of application.

Powerful, Safe and Effective

T-POD[®] provides powerful, fast and safe simultaneous circumferential compression of the pelvic region.

The photos below demonstrate the **T-POD**[®]'s effectiveness in closing and stabilizing the pelvic ring.



Pre-application of **T-POD**[®]



Post-application of **T-POD**[®]

PYNG MEDICAL

7-13511 Crestwood Place
Richmond BC V6V 2E9 Canada
+1-800-349-7964
info@pyng.com www.pyng.com