



Title CSOP 018 Pelvic Trauma

Version No: 1.0

Effective date: 11/02/22

REF PRO400

PROMETHEUS PELVIC SPLINT

Universal Single-Use Pelvic Splint

Description

The Prometheus Pelvic Splint has been designed to fit all patients from the infant to the bariatric adult. It is intended for the emergency management of suspected pelvic fracture(s). It is light, simple to apply, and can be adjusted by a single operator. It is made from soft material and has no internal stitching in order to minimise pressure points and the risk of skin damage.

Instructions for Use

1. Take neoprene band, fold in half with 'fuzzy yellow' surface outside. Place folded band against patient, with the centre of the band in line with the greater trochanter
2. Fold top half of the band down to lie beside the patient's leg
3. Perform a controlled roll to pass the band underneath the patient to the midline
4. Roll the other way to retrieve the folded band
5. Ensure that the centre of the band is still aligned with the greater trochanters
6. Wrap one end of the neoprene band around the patient
7. Attach the blue triangular anchor to the outer surface of the neoprene band
8. Ensure that the centre of the edge of the triangle is directly over the greater trochanter
9. Cut excess neoprene at the level indicated on the triangular anchor. This will allow greater access to the inguinal region
10. Repeat steps 6-9 on the other side
11. Ensure buckle is central and apply tension to the two blue tapes until sufficient force has been applied
12. Secure blue tapes to the neoprene band to maintain desired tension. Record time of splint application.

NSN: 6515-99-174-5924

Rev 1.0 - 01/15



Prometheus Medical Ltd.
The Old Rectory,
Hope under Dinmore,
Herefordshire, HR6 0PW,
United Kingdom
enquiries@prometheusmed.com

