

ANNEX 2: Diagram of CELOX Z-gauze application



1. Ensure someone is applying direct pressure to the wound with a gauze swab. Tear open the CELOX gauze packet and take hold of one end with the other hand.



2. Tightly pack the unfolding CELOX gauze directly to the source of bleeding. Pack remaining wound cavity with CELOX gauze.



3. Apply firm pressure to the wound for THREE minutes. If bleeding persists apply direct pressure for an additional three minutes.



4. Wrap and tie with a bandage/emergency dressing so as to maintain pressure on the wound.